

Charge Plate Holders

Carstens metal or plastic Charge Plate Holders for keeping patient charge plates with the chart can be factory-installed in a variety of positions to fit your needs. Either Charge Plate Holder fits on 1"-4" ringbinders.



Charge Plate Attachment Positions For Ringbinder Charts

Call Customer Service for quote. Factory-installed attachments are not returnable. (See policy on page 107)

Cat. No.	Description	Quantity/Price			
		1-35	36-71	72-239	240+
Charge Plate Options					
1613-01	Plastic Charge Plate Holder	5.81	5.59	5.34	5.13
1611-00	Metal Charge Plate Holder	8.59	8.24	7.89	7.58

Plastic	Metal	
<p>Positions 3C and 4C For Side-Opening Ringbinders</p>	<p>Position 1F and 2F For Side-Opening Ringbinders</p>	<p>Positions 1E and 2E For Top-Opening Ringbinders</p>
<p>Positions 3B and 4B For Top-Opening Ringbinders</p>	<p>Positions 3F and 4F For Side-Opening Ringbinders</p>	<p>Positions 3E and 4E For Top-Opening Ringbinders</p>

Charge Plate Attachment Positions For Springloaded Chartholders

Plastic	Metal		
<p>Position 2A For Springloaded Chartholders</p>	<p>Position 1D and 2D For Springloaded Chartholders</p>		



Desk Charge Plate Storage

Perfect for storing patient charge plates at the nursing station. Holds 15 charge plates each.
Cat. No. 2661-01 Almond
Cat. No. 2661-02 Gray



Charge Plate Storage Rack

Mount on wall or use easel for countertop access. Holds 20 charge plates
Cat. No. 2662-00 White

Cat. No.	Color	Description	Size	Quantity/Price
Charge Plate Storage				
2661-01	Almond	Desk Charge Plate Storage	5½"W x 9¼"L x 4½"H	33.00 each
2661-02	Gray	Desk Charge Plate Storage	5½"W x 9¼"L x 4½"H	33.00 each
2662-00	White	Charge Plate Storage Rack	5¾"W x 14"L x 1"D	93.75 each